



WELLSTOX

long
experience

special
skill

competitive
price

+

+

+

+

+

Lifting Threads

WELLSTOX <

WELLSTOX MONO

- > General Type Thread
- > None Cog
- > Normally use for making volume
- > Easy for control



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E MONO 31-2530-70	31G	25	30	7/0	○	○	○	○
E MONO 30-2530-60	30G	25	30	6/0	○	○	○	○
E MONO 30-3850-60		38	50	6/0	○	○	○	○
E MONO 29-2530-60	29G	25	30	6/0	○	○	○	○
E MONO 29-3850-60		38	50	6/0	○	○	○	○
E MONO 29-5070-60		50	70	6/0	○	○	○	○
E MONO 29-6080-60		60	80	6/0	○	○	○	○
E MONO 27-3850-50	27G	38	50	5/0	○	○	○	○
E MONO 27-5070-50		50	70	5/0	○	○	○	○
E MONO 27-6080-50		60	80	5/0	○	○	○	○
E MONO 27-90150-50		90	150	5/0	○	○	○	○
E MONO 26-3850-50	26G	38	50	5/0	○	○	○	○
E MONO 26-5070-50		50	70	5/0	○	○	○	○
E MONO 26-6080-50		60	80	5/0	○	○	○	○
E MONO 26-90150-50		90	150	5/0	○	○	○	○

Thread detail



WELLSTOX SCREW

- > Twisted thread around needle
- > It can bring the volumizing effect by making collagen



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E SCREW 31-2530-70	31G	25	30	7/0	○	○	○	○
E SCREW 30-2530-60	30G	25	30	6/0	○	○	○	○
E SCREW 30-3850-60		38	50	6/0	○	○	○	○
E SCREW 29-2530-60	29G	25	30	6/0	○	○	○	○
E SCREW 29-3850-60		38	50	6/0	○	○	○	○
E SCREW 29-5070-60		50	70	6/0	○	○	○	○
E SCREW 29-6080-60		60	80	6/0	○	○	○	○
E SCREW 27-3850-50	27G	38	50	5/0	○	○	○	○
E SCREW 27-5070-50		50	70	5/0	○	○	○	○
E SCREW 27-6080-50		60	80	5/0	○	○	○	○
E SCREW 26-6080-50	26G	60	80	5/0	○	○	○	○
E SCREW 26-90150-50		90	150	5/0	○	○	○	○

Thread detail



WELLSTOX MULTI

- > It looks only one filament but made by lots of filaments
- > None Cog
- > Natural volume effect



CODE	N/G	N/L	T/L	USP	PGA
E MULTI PGA 30-2530-60	30G	25	30	6/0	○
E MULTI PGA 30-3850-60		38	50	6/0	○
E MULTI PGA 29-2530-60	29G	25	30	6/0	○
E MULTI PGA 29-3850-60		38	50	6/0	○
E MULTI PGA 29-5070-60		50	70	6/0	○
E MULTI PGA 27-3850-50	27G	38	50	5/0	○
E MULTI PGA 27-5070-50		50	70	5/0	○
E MULTI PGA 27-6080-50		60	80	5/0	○
E MULTI PGA 26-3850-50	26G	38	50	5/0	○
E MULTI PGA 26-5070-50		50	70	5/0	○
E MULTI PGA 26-6080-50		60	80	5/0	○

Thread detail



WELLSTOX TWIN

- > 2 Threads are twisted
- > None Cog
- > Natural volume effect



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TWIN SCREW 27-3850-70	27G	38	50	7/0	○	○	○	
E TWIN SCREW 27-5070-70		50	70	7/0	○	○	○	
E TWIN SCREW 27-6080-70		60	80	7/0	○	○	○	
E TWIN SCREW 26-3850-70	26G	38	50	7/0	○	○	○	
E TWIN SCREW 26-5070-70		50	70	7/0	○	○	○	
E TWIN SCREW 26-6080-70		60	80	7/0	○	○	○	

WELLSTOX G LINE

- > Good for eye area wrinkle and whitening
- > Semipermanent
- > Gold 99.99%



Thread detail




CODE	N/G	N/L	T/L	USP	GOLD
G Line 30-1315-50	30G	13	15	5/0	○
G Line 30-2532 5-50		25	32.5	5/0	○
G Line 30-3850-50		38	50	5/0	○
G Line 29-2532 5-50	29G	25	32.5	5/0	○
G Line 29-3850-50		38	50	5/0	○

WELLSTOX TWIN SCREW

- > Twisted thread (2pcs) are twisted one more outside of needle
- > None Cog
- > Natural volume effect



Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TWIN SCREW 27-3850-70	27G	38	50	7/0	○	○	○	
E TWIN SCREW 27-5070-70		50	70	7/0	○	○	○	
E TWIN SCREW 27-6080-70		60	80	7/0	○	○	○	
E TWIN SCREW 26-3850-70	26G	38	50	7/0	○	○	○	
E TWIN SCREW 26-5070-70		50	70	7/0	○	○	○	
E TWIN SCREW 26-6080-70		60	80	7/0	○	○	○	

WELLSTOX TRIPLE

- > 3 Threads are twisted
- > None Cog
- > Natural volume effect




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TRIPLE 27-3850-70	27G	38	50	7/0	○			
E TRIPLE 27-5070-70		50	70	7/0	○			
E TRIPLE 27-6080-70		60	80	7/0	○			
E TRIPLE 26-3850-70	26G	38	50	7/0	○			
E TRIPLE 26-5070-70		50	70	7/0	○			
E TRIPLE 26-6080-70		60	80	7/0	○			

WELLSTOX BOOM-6

- > Made with 6 pcs of filaments
- > None Cog
- > Natural volume effect




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E BOOM-6 23-3850-60	23G	38	50	6/0 6EA	○	○	○	○
E BOOM-6 23-6080-60		60	80	6/0 6EA	○	○	○	○

WELLSTOX BOOM-10

- > Made with 10 pcs of filaments
- > None Cog
- > Natural volume effect




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E BOOM-10 23-3850-70	23G	38	50	7/0 10EA	○	○	○	○
E BOOM-10 23-6080-70		60	80	7/0 10EA	○	○	○	○

WELLSTOX BOOM-7

- > Made with 7 pcs of filaments
- > None Cog
- > Natural volume effect




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E BOOM-7 23-3850-60	23G	38	50	6/0 7EA	○	○	○	○
E BOOM-7 23-6080-60		60	80	6/0 7EA	○	○	○	○

WELLSTOX NANO

- > Thin but it's exquisite cog thread
- > Use for eye area



Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E N 25-38136-50	25G	38	136	5/0	○			
E N 25-60185-50		60	185	5/0	○			

WELLSTOX CANNULA L

- > Cog thread with L Shape(edge) cannula
- > Painless injection
- > Suitable for the people who can get some trauma during injection




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E CL 23-60115-30	23G	60	115	3/0	○	○	○	
E CL 23-90165-30		90	165	3/0	○	○	○	
E CL 21-60115-20	21G	60	115	2/0	○	○	○	
E CL 21-90165-20		90	165	2/0	○	○	○	
E CL 19-100155-0	19G	100	155	0	○	○	○	
E CL 18-100155-2	18G	100	155	2	○	○	○	

WELLSTOX M1

- > The thread that is made by pressing
- > With W and L type cannula
- > Painless injection
- > Strongest effect
- > Best recommended



Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E M1 CW 19-100170-34	19G	100	170	3&4	○			○
E M1 CL 19-100175-34		100	175	3&4	○			○
E M1 CW 18-100170-34	18G	100	170	3&4	○			○
E M1 CL 18-100175-34		100	175	3&4	○			○

WELLSTOX CANNULA W

- > Cog thread with U Shape(edge) cannula
- > Painless injection
- > Suitable for the people who can get some trauma during injection




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E CW 21-60135-20	21G	60	135	2/0	○	○	○	
E CW 21-90140-20		90	140	2/0	○	○	○	
E CW 19-100150-1	19G	100	150	1	○	○	○	
E CW 18-100150-2	18G	100	150	2	○	○	○	

WELLSTOX M2

- > Made by none heating pressing
- > With W type cannula, painless injection
- > Strongest effect for whole face




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E M2 CW 19-100170-34	19G	100	170	3&4	○			
E M2 CL 19-100175-34		100	175	3&4	○			
E M2 CW 18-100170-34	18G	100	170	3&4	○			
E M2 CL 18-100175-34		100	175	3&4	○			

WELLSTOX MESH COVER

- > Made with mesh type filament
- > None Cog
- > Natural volume effect




Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E MESH COVER CW 23-6025-60-16	23G	60	25	6/0	○		○	
E MESH COVER CW 23-6040-60-16		60	40	6/0	○		○	
E MESH COVER CW 21-6025-60-16	21G	60	25	6/0	○		○	
E MESH COVER CW 21-6040-60-16		60	40	6/0	○		○	

WELLSTOX CLL

- > Both side cog (inside and outside of needle), long thread
- > With cannula L type
- > Good for making face line




Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E CLL 23-60150-30	23G	60	150	3/0	○		○	
E CLL 23-90210-30		90	210	3/0	○		○	
E CLL 21-60150-20	21G	60	150	2/0	○		○	
E CLL 21-90210-20		90	210	2/0	○		○	
E CLL 19-100230-0	19G	100	230	0	○		○	
E CLL 18-100230-2	18G	100	230	2	○		○	

WELLSTOX COG MESH

- > Cog and mesh combined
- > Inside of mesh, there are regeneration tissue with cog -> long effect



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E COG MESH CW 18-100155-2	18G	100	155	2	○			
E COG MESH CL 18-100155-2		100	155	2	○			

WELLSTOX CSL

- > Hold and pull with two pieces of thread at once




Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E CSL 23-38188-30	23G	38	188	3/0	○			
E CSL 23-60240-30		60	240	3/0	○			

WELLSTOX RHINO

- > Specialized thread for elevation of nose tip and bridge




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E R 21-6080-0	21G	60	80	0	○	○	○	
E R 19-3840-2	19G	38	40	2	○	○	○	
E R 19-3850-2		38	50	2	○	○	○	
E R 19-5070-2		50	70	2	○	○	○	
E R 19-6080-2		60	80	2	○	○	○	

WELLSTOX DOUBLE ARM

- > Both side of thread, there are attached 19G or 21G needle (NO need hook for procedure)
- > Our best recommendation thread




Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E LN 21-100300-0	21G	100	300	0	○	○	○	○
E LN 21-120400-0		120	400	0	○	○	○	○
E LN 19-50200-0	19G	50	200	0	○	○	○	○
E LN 19-100300-0		100	300	0	○	○	○	○
E LN 19-120400-1		120	400	1	○	○	○	○
E LN 19-120420-0		120	420	0	○	○	○	○

WELLSTOX E

- > None Cog thread with U Shape(edge) cannula
- > Painless injection
- > Effective for eye dark circles (periorbital melanosis)




Thread detail




CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	HA
E 30-2535-60	30G	25	35	6/0	○	○	○	○
E 30-3850-60		38	50	6/0	○	○	○	○

WELLSTOX DOUBLE ARM M1

- > Same as Double arm but use pressing cog thread (M1)
- > Strong effect, best recommendation



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E LN M1 19-100320-2	19G	100	320	2	○			

WELLSTOX TM1

- > 360° directions rotated pressing cog
- > Hold and pull in 360° directions by pressing cog



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TM1 CW 18-100170-34	18G	100	170	3&4	○			○

WELLSTOX TM3

- > 360° directions rotated pressing cog
- > Hold and pull in 360° directions by pressing cog
- > Same with TM1, the cogs are bigger



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TM3 CW 18-100170-2	18G	100	170	2				○

WELLSTOX TM5

- > 360° directions rotated pressing cog
- > Hold and pull in 360° directions by pressing cog
- > Make stronger support due to larger cog on one side



Thread detail



CODE	N/G	N/L	T/L	USP	PDO	PLLA	PCL	PCLLA
E TM5 CW 18-100170-2	18G	100	170	2				○